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PRESCRIBING INFORMATION

Patients should be counseled that this product does not protect against HIV infection (AIDS) and other sexually transmitted diseases. The BALERINE, IUB™ SCu300B MIDI Intrauterine Device should be placed and removed only by appropriately trained health care professionals.

Device Description

BALERINE, the IUB™ SCu300B MIDI Intrauterine Device (BALERINE) is a sphere-shaped copper intrauterine device (IUD), measuring about 15mm in diameter. A monofilament blue polypropylene double tailed thread is attached to one end of the PET coated frame, each at least 10.5 cm in length, to aid in detection and removal of the device. The BALERINE also contains copper: the total exposed copper surface area is 300mm2. One BALERINE weighs less than one (1) gram. No component of BALERINE or its packaging contains latex. The BALERINE is supplied sterile, pre-loaded in a scaled insertion tube. Also supplied are a solid purple push rod and a sliding flange which is fitted on the insertion tube and aids in gauging the depth of insertion through the cervical canal and into the uterine cavity.

BALERINE - Clinical Pharmacology

The contraceptive effectiveness of the BALERINE is based on the accepted mechanism of action of copper intrauterine devices in which copper has been shown to have spermicidal properties, thus preventing fertilization.

Indications and Usage for BALERINE

The BALERINE is indicated for intrauterine contraception for up to 5 years.

Contraindications

- The BALERINE should not be used in patients younger than 15 years or in patients with the following known or suspected diseases:
1. Pregnancy or suspicion of pregnancy.
 2. Abnormalities of the uterus resulting in distortion of the uterine cavity.
 3. Acute pelvic inflammatory disease or current behavior suggesting a high risk for pelvic inflammatory disease.
 4. Postpartum endometritis or post-abortal endometritis in the past 3 months.
 5. Known or suspected uterine or cervical malignancy.
 6. Genital bleeding of unknown etiology.
 7. Mucopurulent cervicitis - untreated acute cervicitis or vaginitis, including bacterial vaginosis or other lower genital tract infections until infection is controlled.
 8. Wilson’s disease.
 9. Allergy to any component of BALERINE.
 10. A previously placed IUD that has not been removed.
 11. Conditions associated with increased susceptibility to pelvic infection.
 - 12.As generally applied to copper IUDs discretion is advised when considering the use of the BALERINE in women with known anemia, use of anticoagulants, dysmenorrhea or severe menorrhagia.

Warnings and Precautions

The following adverse events are associated with the use of copper intrauterine devices:

1. Intrauterine pregnancy.
2. Ectopic pregnancy.
3. Pelvic infection.
4. Embedment.
5. Perforation.
6. Mal-positioning.
7. Wilson’s disease.
8. Vaginal bleeding.
9. Vasovagal reactions, including fainting during or right after insertion.
10. Expulsion following routine placement and especially after birth or abortion.
11. Magnetic resonance imaging (MRI) personnel must be made aware of the presence of the BALERINE prior to imaging.

12. Medical diathermy.
13. Septic abortion.
14. Sexual partner discomfort stemming from untrimmed removal threads.

Side Effects

Adverse events associated with intrauterine contraception are discussed in WARNINGS and PRECAUTIONS.

INSTRUCTIONS FOR USE

The placement technique for the BALERINE is similar to that used for other IUDs. The health care professional should be familiar with the following instructions. The BALERINE may be placed at any time during the cycle when the appropriately trained health care professional is reasonably certain the patient is not pregnant. However, it is preferable to insert the BALERINE during menstruation. A single BALERINE should be placed in the uterine cavity. The BALERINE should be removed on or before 5 years from the date of insertion.

Before Placement:

1. Make sure that the patient is an appropriate candidate for BALERINE and that she has read the Patient information leaflet.
2. Use of an analgesic before insertion is at the discretion of the patient and the clinician.
3. Establish the size and position of the uterus by pelvic examination.
4. Insert a speculum and cleanse the vagina and cervix with an antiseptic solution.
5. Application of a tenaculum to the cervix is optional for gentle traction of the cervical canal with the uterine cavity.
6. Gently insert a sterile sound to measure the depth of the uterine cavity.
7. The uterus should sound to a depth of 6 to 9 cm except when inserting the BALERINE immediately post-abortion or post-partum. Insertion of BALERINE into a uterine cavity measuring less than 6 cm may increase the incidence of expulsion, bleeding, pain, and perforation. If you encounter cervical stenosis, avoid undue force. Dilators may be helpful in this situation.

How to Place the BALERINE:

- STEP 1** Open the sterile package; Use the flange to mark on the insertion tube the uterine depth that you measured with the sound. Pass the loaded insertion tube through the cervical canal until the gauge is in touch with the cervical external os.
- STEP 2** Insert the push rod into the insertion tube and push forwards in moderation to deploy the BALERINE into the uterine cavity. To ensure the BALERINE is properly positioned avoid insertion tube withdrawal before or during deployment.
- STEP 3** Pull out the rod and then pull out the insertion tube.
- STEP 4** Cut the threads about 2cm out of the cervix.
- STEP 5** **Perform an US examination to make sure the BALERINE is in the center of uterine cavity.** If the BALERINE is not positioned completely within the uterus, remove it and replace it with a new BALERINE. Do not reinsert an expelled or partially expelled BALERINE

CAUTION

Instrumentation of the cervical os may result in vasovagal reactions, including fainting. Have the patient remain supine until she feels well, and have her get up with caution.

Continuing Care

Following placement, examine the patient after her first menses to confirm that the BALERINE is still in place. You should be able to see or feel only the threads. If the BALERINE has been partially or completely expelled, remove it. You can place a new BALERINE if the patient desires and if she is not pregnant. Do not reinsert a used BALERINE.


- Evaluate the patient promptly if she complains of any of the following:
- Abdominal or pelvic pain, cramping, or tenderness; malodorous discharge; bleeding; fever.
 - A missed period.

The length of the visible threads may change with time. However, no action is needed unless you suspect partial expulsion, perforation, or pregnancy.

If you cannot find the threads in the vagina, check that the BALERINE is still in the uterus. The threads can retract into the uterus or break, or the BALERINE may have perforated the uterus or expelled. Radiography or sonography may be required to locate the BALERINE. If there is evidence of partial expulsion, perforation or breakage, remove the BALERINE.

How to Remove the BALERINE

Remove the BALERINE with forceps, pulling gently on the exposed threads. Inspect to assure the integrity of the BALERINE, specifically to the presence of the leading and trailing copper balls (see image below). You may immediately insert a new BALERINE if the patient requests this and has no contraindications.

Embedment or breakage of the BALERINE in the myometrium can make removal difficult. Analgesia, paracervical anesthesia, and cervical dilation may assist in removing an embedded BALERINE. An alligator forceps or other grasping instrument may be helpful. Hysteroscopy may also be helpful.

How is the BALERINE Supplied

The BALERINE is available in cartons of 1 (one) sterile unit. Each BALERINE is packaged in a sterile pouch together with an insertion tube, a slider and a push rod. The BALERINE is supplied sterile. Method of sterilization is ethylene oxide.

Shipping and Storage Conditions

Store the BALERINE packaging in a dry environment at 15°C to 30°C. In these conditions the BALERINE’s shelf life is 3 years. Short term transportation of the BALERINE packaging should be limited to a temperature between -18°C and 55°C.

- SINGLE USE, DO NOT RSTERILIZE.
- NEVER RE-INSERT A USED BALERINE.
- NEVER USE A BALERINE IF THE PACKAGE IS DAMAGED OR OPEN.
- DO NOT USE PAST THE EXPIRY DATE.
- DISPOSE OF USED BALERINE AND ITS COMPONENTS USING BIO-HAZARD DISPOSAL PRACTICES.

- Single use only.
- Ethylene oxide sterilized.
- Do not resterilize.
- Do not use if package is damaged or open.
- Single unit per package.
- Consult instructions for use.
- Caution. There are specific precautions related to the device, refer to IFU.
- Recycle package after use.

BALERINE is a trademark of OCON Medical Ltd. www.oconmed.com	OCON Medical Ltd. 15 Hashdera Hamerkazit, POB 552, Modiin 7171801, Israel Tel. (972) 72 21 50 105
--------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------

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